

Product datasheet for TA810522

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

KIAA1970 (EARS2) Mouse Monoclonal Antibody [Clone ID: OTI2C1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2C1
Applications: IHC, WB

Recommended Dilution: WB 1:500~2000, IHC 1:500

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombiant protein of human EARS2 (NP_001077083) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 58.5 kDa

Gene Name: glutamyl-tRNA synthetase 2, mitochondrial

Database Link: NP 001077083

Entrez Gene 67417 MouseEntrez Gene 361641 RatEntrez Gene 124454 Human

Q5JPH6

Synonyms: COXPD12; gluRS; MSE1


Protein Families: Druggable Genome

Protein Pathways: Aminoacyl-tRNA biosynthesis, Metabolic pathways, Porphyrin and chlorophyll metabolism






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY EARS2 (Cat# [RC221871], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-EARS2 (Cat# TA810522)(1:2000). Positive lysates [LY421237] (100ug) and [LC421237] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-EARS2 monoclonal antibody (1:500).






Immunohistochemical staining of paraffinembedded Human breast tissue within the normal limits using anti-EARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA810522) (1:500)




Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-EARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA810522) (1:500)




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-EARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA810522) (1:500)





Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-EARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA810522) (1:500)



Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-EARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA810522) (1:500)