

Product datasheet for **TA810580S**

RNF149 Mouse Monoclonal Antibody [Clone ID: OTI8D7]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI8D7
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 250-383 of human RNF149 (NP_775918) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	43 kDa
Gene Name:	ring finger protein 149
Database Link:	NP_775918 Entrez Gene 284996 Human Q8NC42
Synonyms:	DNAPTP2
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)



Product images:


Figure A, Western blot analysis of overexpressed lysates (15ug per lane) from HEK293T cells transfected with empty plasmid ([PS100001], NC), human RNF149 plasmid ([RC207979], hRNF149), mouse RNF149 plasmid ([MR218021], m RNF149) using anti- RNF149 antibody [TA810580](1:500). Figure B, Western blot analysis of the same samples as figure A with anti-DDK antibody ([TA180144], 1:1000)



Western blot analysis of extracts (35ug) from 5 cell lines by using anti-RNF149 monoclonal antibody. (1:2000) (1:500)