


## Product datasheet for **TA810658M**

### NEUROD4 Mouse Monoclonal Antibody [Clone ID: OTI1B7]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1B7
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human NEUROD4 (NP_067014) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	neuronal differentiation 4
Database Link:	<a href="#">NP_067014</a> <a href="#">Entrez Gene 11923 Mouse</a> <a href="#">Entrez Gene 288821 Rat</a> <a href="#">Entrez Gene 58158 Human</a> <a href="#">Q9HD90</a>
Background:	Probably acts as a transcriptional activator. Mediates neuronal differentiation. Required for the regulation of amacrine cell fate specification in the retina. [UniProtKB/Swiss-Prot Function]
Synonyms:	ATH-3; ATH3; Atoh3; bHLHa4; MATH-3; MATH3
Protein Families:	Druggable Genome, Transcription Factors


[View online »](#)

**Product images:**


HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY NEUROD4 (Cat# [RC222464], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NEUROD4 (Cat# [TA810658])(1:2000). Positive lysates [LY412030] (100ug) and [LC412030] (20ug) can be purchased separately from OriGene.