


## Product datasheet for **TA810669M**

### **ZNF175 Mouse Monoclonal Antibody [Clone ID: OTI1F2]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI1F2
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ZNF175 (NP_009078) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	zinc finger protein 175
Database Link:	<a href="#">NP_009078</a> <a href="#">Entrez Gene 7728 Human</a> <a href="#">Q9Y473</a>
Synonyms:	OTK18
Protein Families:	Transcription Factors

[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ZNF175 (Cat# [RC208740], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF175 (Cat# [TA810669])(1:2000). Positive lysates [LY416153] (100ug) and [LC416153] (20ug) can be purchased separately from OriGene.