

Product datasheet for **TA810684AM**


GRHL3 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI5F3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI5F3
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human GRHL3 (NP_937817) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	56.9 kDa
Gene Name:	grainyhead like transcription factor 3
Database Link:	NP_937817 Entrez Gene 230824 Mouse Entrez Gene 298555 Rat Entrez Gene 57822 Human Q8TE85
Synonyms:	SOM; TFCP2L4; VWS2



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY GRHL3 (Cat# [RC207857], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GRHL3 (Cat# [TA810684])(1:2000). Positive lysates [LY404981] (100ug) and [LC404981] (20ug) can be purchased separately from OriGene.