

## Product datasheet for **TA810707S**


### **RPB11 (POLR2J) Mouse Monoclonal Antibody [Clone ID: OTI1H4]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OTI1H4
<b>Applications:</b>	WB
<b>Recommended Dilution:</b>	WB 1:2000
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Full length human recombinant protein of human POLR2J (NP_006225) produced in HEK293T cell.
<b>Formulation:</b>	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
<b>Concentration:</b>	1 mg/ml
<b>Purification:</b>	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Predicted Protein Size:</b>	13.1 kDa
<b>Gene Name:</b>	polymerase (RNA) II subunit J
<b>Database Link:</b>	<a href="#">NP_006225</a> <a href="#">Entrez Gene 20022 Mouse</a> <a href="#">Entrez Gene 288588 Rat</a> <a href="#">Entrez Gene 5439 Human</a> <a href="#">P52435</a>
<b>Synonyms:</b>	hRPB14; POLR2J1; RPB11; RPB11A; RPB11m
<b>Protein Families:</b>	Transcription Factors
<b>Protein Pathways:</b>	Huntington's disease, Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY POLR2J (Cat# [RC204819], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-POLR2J (Cat# [TA810707]). Positive lysates [LY401876] (100ug) and [LC401876] (20ug) can be purchased separately from OriGene.