

Product datasheet for **TA810772BM**

APOBEC4 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1F7]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1F7
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human APOBEC4 (NP_982279) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	apolipoprotein B mRNA editing enzyme catalytic polypeptide like 4
Database Link:	NP_982279 Entrez Gene 71281 Mouse Entrez Gene 403314 Human Q8WW27
Synonyms:	C1orf169




[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY APOBEC4 ([RC205405], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-APOBEC4 (1:500).



Western blot analysis of extracts (35ug) from HT29 and mouse testis by using anti-APOBEC4 monoclonal antibody (1:500).