

## Product datasheet for **TA810772S**

### **APOBEC4 Mouse Monoclonal Antibody [Clone ID: OT11F7]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OT11F7
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human APOBEC4 (NP_982279) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	apolipoprotein B mRNA editing enzyme catalytic polypeptide like 4
Database Link:	<a href="#">NP_982279</a> <a href="#">Entrez Gene 71281 Mouse</a> <a href="#">Entrez Gene 403314 Human</a> <a href="#">Q8WW27</a>
Synonyms:	C1orf169




[View online »](#)

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY APOBEC4 ([RC205405], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-APOBEC4 (1:500).



Western blot analysis of extracts (35ug) from HT29 and mouse testis by using anti-APOBEC4 monoclonal antibody (1:500).