

## Product datasheet for **TA810800S**


### TAZ Mouse Monoclonal Antibody [Clone ID: OTI5A4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI5A4
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 45-248 of human TAZ (NP_851830) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	tafazzin
Database Link:	<a href="#">NP_851830</a> <a href="#">Entrez Gene 66826 Mouse</a> <a href="#">Entrez Gene 363521 Rat</a> <a href="#">Entrez Gene 6901 Human</a> <a href="#">Q16635</a>
Synonyms:	BTHS; CMD3A; EFE; EFE2; G4.5; LVNCX; TAZ; Taz1
Protein Families:	ES Cell Differentiation/IPS, Transmembrane



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TAZ ([RC219401], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TAZ (1:500).