

## Product datasheet for **TA810842BM**


### DPPA5 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI2F6]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2F6
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human DPPA5 (NP_001020461) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	13.3 kDa
Gene Name:	developmental pluripotency associated 5
Database Link:	<a href="#">NP_001020461</a> <a href="#">Entrez Gene 340168 Human</a> <a href="#">A6NC42</a>
Synonyms:	ESG1



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY DPPA5 (Cat# [RC222509], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DPPA5 (Cat# [TA810842])(1:2000). Positive lysates [LY400405] (100ug) and [LC400405] (20ug) can be purchased separately from OriGene.