

Product datasheet for **TA810898S**

Heat Shock Factor 2 Binding Protein (HSF2BP) Mouse Monoclonal Antibody [Clone ID: OTI2A1]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2A1
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-181 of human HSF2BP (NP_008962) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	37.5 kDa
Gene Name:	heat shock transcription factor 2 binding protein
Database Link:	NP_008962 Entrez Gene 74377 Mouse Entrez Gene 499413 Rat Entrez Gene 11077 Human O75031
Synonyms:	heat shock factor 2 binding protein; heat shock transcription factor 2 binding protein
Protein Families:	Transcription Factors



[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY HSF2BP (Cat# [RC201580], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HSF2BP (Cat# [TA810898])(1:500). Positive lysates [LY416250] (100ug) and [LC416250] (20ug) can be purchased separately from OriGene.