

## Product datasheet for **TA810971BM**


### USP25 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI6H4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI6H4
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human USP25 (NP_037528) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	122 kDa
Gene Name:	ubiquitin specific peptidase 25
Database Link:	<a href="#">NP_037528</a> <a href="#">Entrez Gene 30940 Mouse</a> <a href="#">Entrez Gene 29761 Human</a> <a href="#">Q9UHP3</a>
Synonyms:	USP21
Protein Families:	Protease



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY USP25 (Cat# [RC209755], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-USP25 (Cat# [TA810971])(1:2000). Positive lysates [LY415616] (100ug) and [LC415616] (20ug) can be purchased separately from OriGene.