

## Product datasheet for **TA810979AM**

### **PX19 (PRELID1) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI4A4]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI4A4
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 36-219 of human PRELID1 (NP_037369) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	25 kDa
Gene Name:	PRELI domain containing 1
Database Link:	<a href="#">NP_037369</a> <a href="#">Entrez Gene 27166 Human</a> <a href="#">Q9Y255</a>
Synonyms:	CGI-106; PRELI; PX19; SBB112



[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PRELID1 ([RC202303], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PRELID1 (1:2000). Positive lysates [LY415724] (100ug) and [LC415724] (20ug) can be purchased separately from OriGene.