

## Product datasheet for **TA810979S**


### PX19 (PRELID1) Mouse Monoclonal Antibody [Clone ID: OTI4A4]

#### Product data:

|                         |                                                                                                                         |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Type:           | Primary Antibodies                                                                                                      |
| Clone Name:             | OTI4A4                                                                                                                  |
| Applications:           | WB                                                                                                                      |
| Recommended Dilution:   | WB 1:2000                                                                                                               |
| Reactivity:             | Human, Mouse, Rat                                                                                                       |
| Host:                   | Mouse                                                                                                                   |
| Isotype:                | IgG2b                                                                                                                   |
| Clonality:              | Monoclonal                                                                                                              |
| Immunogen:              | Human recombinant protein fragment corresponding to amino acids 36-219 of human PRELID1 (NP_037369) produced in E.coli. |
| Formulation:            | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                                                    |
| Concentration:          | 1 mg/ml                                                                                                                 |
| Purification:           | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)               |
| Conjugation:            | Unconjugated                                                                                                            |
| Storage:                | Store at -20°C as received.                                                                                             |
| Stability:              | Stable for 12 months from date of receipt.                                                                              |
| Predicted Protein Size: | 25 kDa                                                                                                                  |
| Gene Name:              | PRELI domain containing 1                                                                                               |
| Database Link:          | <a href="#">NP_037369</a><br><a href="#">Entrez Gene 27166 Human</a><br><a href="#">Q9Y255</a>                          |
| Synonyms:               | CGI-106; PRELI; PX19; SBB112                                                                                            |



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PRELID1 ([RC202303], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PRELID1 (1:2000). Positive lysates [LY415724] (100ug) and [LC415724] (20ug) can be purchased separately from OriGene.