

## Product datasheet for **TA810996M**

### ADCK3 (COQ8A) Mouse Monoclonal Antibody [Clone ID: OTI2A7]

#### Product data:


|                         |   |
|-------------------------|---|
| Product Type:           | Primary Antibodies  |
| Clone Name:             | OTI2A7  |
| Applications:           | WB  |
| Recommended Dilution:   | WB 1:2000   |
| Reactivity:             | Human   |
| Host:                   | Mouse   |
| Isotype:                | IgG1  |
| Clonality:              | Monoclonal  |
| Immunogen:              | Human recombinant protein fragment corresponding to amino acids 57-328 of human ADCK3 (NP_064632) produced in E.coli.   |
| Formulation:            | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.  |
| Concentration:          | 1 mg/ml   |
| Purification:           | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)   |
| Conjugation:            | Unconjugated  |
| Storage:                | Store at -20°C as received.   |
| Stability:              | Stable for 12 months from date of receipt.  |
| Predicted Protein Size: | 71.8 kDa  |
| Gene Name:              | aarF domain containing kinase 3   |
| Database Link:          | <a href="#">NP_064632</a><br><a href="#">Entrez Gene 56997 Human</a><br><a href="#">Q8NI60</a>  |
| Background:             | This gene encodes a mitochondrial protein similar to yeast ABC1, which functions in an electron-transferring membrane protein complex in the respiratory chain. It is not related to the family of ABC transporter proteins. Expression of this gene is induced by the tumor suppressor p53 and in response to DNA damage, and inhibiting its expression partially suppresses p53-induced apoptosis. Alternatively spliced transcript variants have been found; however, their full-length nature has not been determined. [provided by RefSeq, Jul 2008] |


[View online »](#)

**Synonyms:** ARCA2; CABC1; COQ8; COQ10D4; SCAR9

**Protein Families:** Protein Kinase

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ADCK3 (Cat# [RC201457], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ADCK3 (Cat# [TA810996]). Positive lysates [LY402769] (100ug) and [LC402769] (20ug) can be purchased separately from OriGene.