

Product datasheet for TA811064

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

p63 (TP63) Mouse Monoclonal Antibody [Clone ID: OTI2A7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2A7
Applications: WB

Recommended Dilution: WB 1:500

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Synthetic peptide corresponding to the N-terminus of Human P40 (ΔNp63) (NP_001108452)

conjugated to KLH. The exact sequence is proprietary.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 44.5 kDa

Gene Name: tumor protein p63

Database Link: NP 001108454

Entrez Gene 246334 RatEntrez Gene 8626 Human

Q9H3D4





Background:

This gene encodes a member of the p53 family of transcription factors. An animal model, p63 -/- mice, has been useful in defining the role this protein plays in the development and maintenance of stratified epithelial tissues. p63 -/- mice have several developmental defects which include the lack of limbs and other tissues, such as teeth and mammary glands, which develop as a result of interactions between mesenchyme and epithelium. Mutations in this gene are associated with ectodermal dysplasia, and cleft lip/palate syndrome 3 (EEC3); splithand/foot malformation 4 (SHFM4); ankyloblepharon-ectodermal defects-cleft lip/palate; ADULT syndrome (acro-dermato-ungual-lacrimal-tooth); limb-mammary syndrome; Rap-Hodgkin syndrome (RHS); and orofacial cleft 8. Both alternative splicing and the use of alternative promoters results in multiple transcript variants encoding different proteins. Many transcripts encoding different proteins have been reported but the biological validity and the full-length nature of these variants have not been determined. [provided by RefSeq, Jul 2008]


Synonyms:

AIS; B(p51A); B(p51B); EEC3; KET; LMS; NBP; OFC8; p40; p51; p53CP; p63; p73H; p73L; RHS;

Protein Families:

Druggable Genome, Transcription Factors

Product images:



Western blot analysis of lysates from HEK293T cells transfected with P40 (Δ Np63) (Lane 2, Cat# [RC225987]), TP63 overexpression plasmid (Lane 3, Cat# [RC208013]) or empty vector plasmid (Lane 1) using anti-P40 (Δ Np63) antibody (Cat# TA811064). (1:1000)