

## Product datasheet for **TA811079**


### **RBJ (DNAJC27) Mouse Monoclonal Antibody [Clone ID: OTI7A4]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI7A4
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human DNAJC27 (NP_057628) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	30.7 kDa
Gene Name:	Dnaj heat shock protein family (Hsp40) member C27
Database Link:	<a href="#">NP_057628</a> <a href="#">Entrez Gene 217378 Mouse</a> <a href="#">Entrez Gene 298859 Rat</a> <a href="#">Entrez Gene 51277 Human</a> <a href="#">Q9NZQ0</a>
Synonyms:	RabJ5; RBJ



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DNAJC27 ([RC205500], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DNAJC27 (1:2000). Positive lysates [LY413917] (100ug) and [LC413917] (20ug) can be purchased separately from OriGene.