


Product datasheet for **TA811080M**

RBJ (DNAJC27) Mouse Monoclonal Antibody [Clone ID: OTI7A5]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI7A5
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human DNAJC27 (NP_057628) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	30.7 kDa
Gene Name:	DnaJ heat shock protein family (Hsp40) member C27
Database Link:	NP_057628 Entrez Gene 217378 Mouse Entrez Gene 298859 Rat Entrez Gene 51277 Human Q9NZQ0
Synonyms:	RabJS; RBJ


[View online »](#)

Product images:


HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DNAJC27 ([RC205500], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DNAJC27 (1:2000). Positive lysates [LY413917] (100ug) and [LC413917] (20ug) can be purchased separately from OriGene.