

## Product datasheet for **TA811167BM**


### SOCS4 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI5H9]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI5H9
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human SOCS4 (NP_955453) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	suppressor of cytokine signaling 4
Database Link:	<a href="#">NP_955453</a> <a href="#">Entrez Gene 122809 Human</a> <a href="#">Q8WXH5</a>
Synonyms:	SOCS7
Protein Families:	Druggable Genome
Protein Pathways:	Insulin signaling pathway, Jak-STAT signaling pathway, Type II diabetes mellitus



[View online »](#)

**Product images:**

Western blot analysis of extracts (35ug) from HELA cell line by using anti-SOCS4 monoclonal antibody (1:500).