

## Product datasheet for **TA811193M**

### AIF (AIFM1) Mouse Monoclonal Antibody [Clone ID: OTI5D2]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI5D2
Applications:	IHC, WB
Recommended Dilution:	WB 1:500~2000, IHC 1:500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 134-483 of human AIFM1 (NP_001124318) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	apoptosis inducing factor, mitochondria associated 1
Database Link:	<a href="#">NP_001124318</a> <a href="#">Entrez Gene 26926 Mouse</a> <a href="#">Entrez Gene 83533 Rat</a> <a href="#">Entrez Gene 9131 Human</a> <a href="#">O95831</a>
Synonyms:	AIF; CMT2D; CMTX4; COWCK; COXPD6; NADMR; NAMSD; PDCD8
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Apoptosis




[View online »](#)


## Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AIFM1 ([RC225349], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AIFM1 (1:2000).



Western blot analysis of extracts (35ug) from 2 different cell lines and mouse liver tissue lysate by using anti-AIFM1 monoclonal antibody (1:500).



Equivalent amounts of cell lysates (10 ug per lane) of wild-type 293T cells (WT, Cat# LC810293T) and AIFM1-Knockout 293T cells (KO, Cat# [LC811853]) were separated by SDS-PAGE and immunoblotted with anti-AIFM1 monoclonal antibody [TA811193], (1:500). Then the blotted membrane was stripped and reprobed with anti-β-actin antibody ([TA811000]) as a loading control.



Immunohistochemical staining of paraffin-embedded Human adult heart tissue within the normal limits using anti-AIFM1 mouse monoclonal antibody. Heat-induced epitope retrieval by EDTA solution buffer pH 8.0 at 120°C for 3 min.