

Product datasheet for TA811208BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

AKAP12 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1B6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1B6
Applications: IHC, WB

Recommended Dilution: WB 1:500~2000, IHC 1:500

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 266-557 of human

AKAP12 (NP_005091) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 191.3 kDa

Gene Name: A-kinase anchoring protein 12

Database Link: NP 005091

Entrez Gene 9590 Human

Q02952


Synonyms: AKAP250; SSeCKS

Protein Families: Druggable Genome






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY AKAP12 (Cat# [RC221620], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AKAP12 (Cat# [TA811208])(1:2000). Positive lysates [LY417515] (100ug) and [LC417515] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 2 different cell lines by using anti-AKAP12 monoclonal antibody (1:500).






Immunohistochemical staining of paraffinembedded Human adult heart tissue within the normal limits using anti-AKAP12 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA811208]) (1:500)




Immunohistochemical staining of paraffinembedded Human embryonic cerebellum within the normal limits using anti-AKAP12 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA811208]) (1:500)



Immunohistochemical staining of paraffinembedded Human spleen tissue within the normal limits using anti-AKAP12 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA811208]) (1:500)





Immunohistochemical staining of paraffinembedded Human testicle tissue within the normal limits using anti-AKAP12 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA811208]) (1:500)