

## Product datasheet for **TA811233AM**

### **E74 like factor 1 (ELF1) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OT11H1]**


#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OT11H1
<b>Applications:</b>	WB
<b>Recommended Dilution:</b>	WB 1:500~2000
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Full length human recombinant protein of human ELF1 (NP_758961) produced in E.coli.
<b>Formulation:</b>	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
<b>Concentration:</b>	0.5 mg/ml
<b>Purification:</b>	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
<b>Conjugation:</b>	Biotin
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Gene Name:</b>	E74 like ETS transcription factor 1
<b>Database Link:</b>	<a href="#">NP_758961</a> <a href="#">Entrez Gene 13709 Mouse</a> <a href="#">Entrez Gene 85424 Rat</a> <a href="#">Entrez Gene 1997 Human P32519</a>
<b>Synonyms:</b>	E74-like factor 1 (ets domain transcription factor); OTTHUMP00000018310
<b>Protein Families:</b>	Transcription Factors



[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ELF1 ([RC207140], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ELF1 (1:2000).