

## Product datasheet for **TA811242AM**


### Sec8 (EXOC4) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1E9]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1E9
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human EXOC4 (NP_001032203) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	53.7 kDa
Gene Name:	exocyst complex component 4
Database Link:	<a href="#">NP_001032203</a> <a href="#">Entrez Gene 20336 Mouse</a> <a href="#">Entrez Gene 116654 Rat</a> <a href="#">Entrez Gene 60412 Human</a> <a href="#">Q96A65</a>
Synonyms:	SEC8; SEC8L1; Sec8p
Protein Pathways:	Tight junction



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY EXOC4 ([RC221799], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-EXOC4 (1:500). Positive lysates [LY422147] (100ug) and [LC422147] (20ug) can be purchased separately from OriGene.