

Product datasheet for **TA811266S**

ZNF19 Mouse Monoclonal Antibody [Clone ID: OTI10C3]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI10C3
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ZNF19 (NP_008892) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	zinc finger protein 19
Database Link:	NP_008892 Entrez Gene 7567 Human P17023
Synonyms:	KOX12
Protein Families:	Transcription Factors




[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ZNF19 (Cat# [RC208208], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF19 (Cat# [TA811266])(1:2000).



Western blot analysis of extracts (35ug) from A549 cell line by using anti-ZNF19 monoclonal antibody (1:500).