

## Product datasheet for **TA811299**


### **C3orf62 Mouse Monoclonal Antibody [Clone ID: OTI9F5]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI9F5
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human C3ORF62 (NP_940964) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	30 kDa
Gene Name:	chromosome 3 open reading frame 62
Database Link:	<a href="#">NP_940964</a> <a href="#">Entrez Gene 375341 Human</a> <a href="#">Q6ZUJ4</a>
Synonyms:	FLJ43654; MGC23381; MGC61663; MGC62079



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY C3orf62 (Cat# [RC204956], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-C3orf62 (Cat# TA811299)(1:2000). Positive lysates [LY404856] (100ug) and [LC404856] (20ug) can be purchased separately from OriGene.