

Product datasheet for **TA811359BM**

RNA5-8SN2 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI10D7]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI10D7
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant full length E6&E7 protein of HPV type 18 produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	RNA, 5.8S ribosomal N2



[View online »](#)

Product images:



Western blot analysis of extracts (35ug) from three different cell lysates by using anti-HPV18 E7 monoclonal antibody (1:500).