

## Product datasheet for **TA811372M**


### **TAF1 48 (TAF1A) Mouse Monoclonal Antibody [Clone ID: OTI7A11]**

#### **Product data:**


Product Type:	Primary Antibodies
Clone Name:	OTI7A11
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human TAF1A (NP_647603) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	39.4 kDa
Gene Name:	TATA-box binding protein associated factor, RNA polymerase I subunit A
Database Link:	<a href="#">NP_647603</a> <a href="#">Entrez Gene 9015 Human</a> <a href="#">Q15573</a>
Synonyms:	MGC:17061; RAFI48; SL1; TAFI48
Protein Families:	Transcription Factors


[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY TAF1A (Cat# [RC203755], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TAF1A (Cat# [TA811372])(1:500). Positive lysates [LY408297] (100ug) and [LC408297] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from HELA cell line by using anti-TAF1A monoclonal antibody (1:500).