

Product datasheet for **TA811428BM**

FAM54A (MTFR2) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI6C8]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI6C8
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FAM54A (NP_001092756) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	43.2 kDa
Gene Name:	mitochondrial fission regulator 2
Database Link:	NP_001092756 Entrez Gene 113115 Human Q6P444
Synonyms:	DUFD1; FAM54A



[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FAM54A ([RC213401], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FAM54A (1:2000). Positive lysates [LY426058] (100ug) and [LC426058] (20ug) can be purchased separately from OriGene.