

## Product datasheet for **TA811490AM**


### ZNF181 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI9G2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI9G2
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ZNF181 (NP_001025168) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	zinc finger protein 181
Database Link:	<a href="#">NP_001025168</a> <a href="#">Entrez Gene 339318 Human</a> <a href="#">Q2M3W8</a>
Synonyms:	HHZ181



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ZNF181 (Cat# [RC207393], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF181 (Cat# [TA811490])(1:2000).