

## Product datasheet for TA811528AM

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### CD14 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI12C8]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI12C8

Applications: FC, IHC, WB

Recommended Dilution: WB 1:500, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 20-367 of human CD14

(NP\_000582) produced in E.coli.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 40.08 kDa

**Gene Name:** CD14 molecule

Database Link: NP 000582

Entrez Gene 929 Human

P08571

**Background:** The protein encoded by this gene is a surface antigen that is preferentially expressed on

monocytes/macrophages. It cooperates with other proteins to mediate the innate immune response to bacterial lipopolysaccharide. Alternative splicing results in multiple transcript

variants encoding the same protein. [provided by RefSeq, Mar 2010]

Synonyms: CD14 antigen; CD14 molecule; monocyte differentiation antigen CD14






Protein Families: Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS,

Transmembrane


Protein Pathways: Hematopoietic cell lineage, MAPK signaling pathway, Pathogenic Escherichia coli infection,

Regulation of actin cytoskeleton, Toll-like receptor signaling pathway

# **Product images:**




Western blot analysis of extracts (35ug) from PBMC by using anti-CD14 monoclonal antibody (1:500).




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-CD14 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA811528]) (1:2000)






Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-CD14 mouse monoclonal antibody. (Heatinduced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA811528]) (1:2000)




Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-CD14 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA811528]) (1:2000)




Immunohistochemical staining of paraffinembedded Human Gastric Carcinoma using anti-CD14 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA811528]) (1:2000)





Immunohistochemical staining of paraffinembedded Human spleen tissue within the normal limits using anti-CD14 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA811528]) (1:2000)



HEK293T cells transfected with either [RC203819] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-CD14 antibody ([TA811528]), and then analyzed by flow cytometry (1:100).