

## Product datasheet for **TA811551**


### **FOLR3 Mouse Monoclonal Antibody [Clone ID: OTI6C1]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI6C1
Applications:	IHC
Recommended Dilution:	IHC 1:150
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FOLR3 (NP_000795) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	folate receptor 3 (gamma)
Database Link:	<a href="#">NP_000795</a> <a href="#">Entrez Gene 2352 Human P41439</a>
Synonyms:	FR-G; FR-gamma; gamma-hFR
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

**Product images:**

Immunohistochemical staining of paraffin-embedded Human spleen tissue within the normal limits using anti-FOLR3 mouse monoclonal antibody. Heat-induced epitope retrieval by EDTA solution buffer pH 8.0 at 120°C for 3 min.