

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA811601AM

CRMP2 (DPYSL2) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI5D7]

Product data:

| Product Type: | Primary Antibodies |
|-----------------------|--|
| Clone Name: | OTI5D7 |
| Applications: | WB |
| Recommended Dilution: | WB 1:500~2000 |
| Reactivity: | Human, Rat, Mouse |
| Host: | Mouse |
| lsotype: | lgG2b |
| Clonality: | Monoclonal |
| Immunogen: | Human recombinant protein fragment corresponding to amino acids 13-490 of human DPYSL2 (NP_001377) produced in E.coli. |
| Formulation: | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide. |
| Concentration: | 0.5 mg/ml |
| Purification: | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Conjugation: | Biotin |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Gene Name: | dihydropyrimidinase like 2 |
| Database Link: | <u>NP_001377</u> <u>Entrez Gene 12934 MouseEntrez Gene 25416 RatEntrez Gene 1808 Human</u> <u>Q16555</u> |



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


CRMP2 (DPYSL2) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI5D7] – TA811601AM

Background:This gene encodes a member of the collapsin response mediator protein family. Collapsin
response mediator proteins form homo- and hetero-tetramers and facilitate neuron
guidance, growth and polarity. The encoded protein promotes microtubule assembly and is
required for Sema3A-mediated growth cone collapse, and also plays a role in synaptic
signaling through interactions with calcium channels. This gene has been implicated in
multiple neurological disorders, and hyperphosphorylation of the encoded protein may play
a key role in the development of Alzheimer's disease. Alternatively spliced transcript variants
encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Sep 2011]


| Synonyms: | CRMP-2; CRMP2; DHPRP2; DRP-2; DRP2; N2A3; ULIP-2; ULIP2 |
|-------------------|---|
| Protein Families: | Druggable Genome |

Protein Pathways: Axon guidance

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY DPYSL2 (Cat# [RC209080], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DPYSL2 (1:2000) (Cat# [TA811601]). Positive lysates [LY419992] (100ug) and [LC419992] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-DPYSL2 monoclonal antibody (1:500).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US