

## Product datasheet for **TA811667BM**

### ZNF572 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI2F6]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2F6
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ZNF572 (NP_689625) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	zinc finger protein 572
Database Link:	<a href="#">NP_689625</a> <a href="#">Entrez Gene 137209 Human</a> <a href="#">Q7Z3I7</a>
Synonyms:	FLJ38002




[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ZNF572 (Cat# [RC224676], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF572 (Cat# [TA811667])(1:500). Positive lysates [LY407580] (100ug) and [LC407580] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from HepG2 cell line by using anti-ZNF572 monoclonal antibody (1:500).