

## Product datasheet for **TA811699M**

### **CAPON (NOS1AP) Mouse Monoclonal Antibody [Clone ID: OTI10F8]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI10F8
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-196 of human NOS1AP (NP_055512) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	56 kDa
Gene Name:	nitric oxide synthase 1 adaptor protein
Database Link:	<a href="#">NP_055512</a> <a href="#">Entrez Gene 9722 Human</a> <a href="#">O75052</a>




[View online »](#)

**Background:**

This gene encodes a cytosolic protein that binds to the signaling molecule, neuronal nitric oxide synthase (nNOS). This protein has a C-terminal PDZ-binding domain that mediates interactions with nNOS and an N-terminal phosphotyrosine binding (PTB) domain that binds to the small monomeric G protein, Dexras1. Studies of the related mouse and rat proteins have shown that this protein functions as an adapter protein linking nNOS to specific targets, such as Dexras1 and the synapsins. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2009]

**Synonyms:**

6330408P19Rik; CAPON; NPHS22

**Product images:**


HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY NOS1AP (Cat# [RC216122], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NOS1AP (Cat# [TA811699])(1:2000). Positive lysates [LY402366] (100ug) and [LC402366] (20ug) can be purchased separately from OriGene.