

## Product datasheet for **TA811801**

### Ovary specific acidic protein (MGARP) Mouse Monoclonal Antibody [Clone ID: OTI3F3]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3F3
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human MGARP (NP_116012) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	25.2 kDa
Gene Name:	mitochondria localized glutamic acid rich protein
Database Link:	<a href="#">NP_116012</a> <a href="#">Entrez Gene 84709 Human</a> <a href="#">Q8TDB4</a>
Synonyms:	C4orf49; CESP-1; HUMMR; OSAP
Protein Families:	Transmembrane



[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MGARP ([RC216908], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MGARP (1:2000).