

Product datasheet for **TA811916BM**

KAL1 (ANOS1) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4B8]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI4B8
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 129-373 of human KAL1(NP_000207) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	76.11 kDa
Gene Name:	anosmin 1
Database Link:	NP_000207 Entrez Gene 3730 Human P23352
Background:	Mutations in this gene cause the X-linked Kallmann syndrome. The encoded protein is similar in sequence to proteins known to function in neural cell adhesion and axonal migration. In addition, this cell surface protein is N-glycosylated and may have anti-protease activity. [provided by RefSeq, Jul 2008]
Synonyms:	ADMLX; HH1; HHA; KAL; KAL1; KALIG-1; KMS; WFDC19




[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY KAL1 (Cat# [RC212305], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-KAL1 (Cat# [TA811916]). Positive lysates [LY400084] (100ug) and [LC400084] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from MDA-MB453 cell line by using anti-KAL1 monoclonal antibody (1:500).