

## Product datasheet for **TA812059M**

### EPO Mouse Monoclonal Antibody [Clone ID: OTI5D8]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI5D8
Applications:	WB
Recommended Dilution:	WB 1:150
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 28-193 of human EPO (NP_000790) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	21.31 kDa
Gene Name:	erythropoietin
Database Link:	<a href="#">NP_000790</a> <a href="#">Entrez Gene 2056 Human</a> <a href="#">P01588</a>
Background:	This gene is a member of the EPO/TPO family and encodes a secreted, glycosylated cytokine composed of four alpha helical bundles. The protein is found in the plasma and regulates red cell production by promoting erythroid differentiation and initiating hemoglobin synthesis. This protein also has neuroprotective activity against a variety of potential brain injuries and antiapoptotic functions in several tissue types. [provided by RefSeq, Jul 2008]



[View online »](#)

<b>Synonyms:</b>	EP; MVCD2
<b>Protein Families:</b>	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein, Transmembrane
<b>Protein Pathways:</b>	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway

**Product images:**


HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY EPO (Cat# [RC210775], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-EPO (Cat# [TA812059]). Positive lysates [LY400277] (100ug) and [LC400277] (20ug) can be purchased separately from OriGene.