

# **Product datasheet for TA812166BM**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## ATP6V0D2 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI2D6]

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI2D6
Applications: IHC, WB

Recommended Dilution: WB 1:500~2000, IHC 1:2000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ATP6V0D2 (NP\_689778) produced in E.coli.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

Predicted Protein Size: 40.2 kDa

**Gene Name:** ATPase H+ transporting V0 subunit d2

Database Link: NP 689778

Entrez Gene 245972 Human

Q8N8Y2

Synonyms: ATP6D2; VMA6


**Protein Pathways:** Epithelial cell signaling in Helicobacter pylori infection, Lysosome, Metabolic pathways,

Oxidative phosphorylation, Vibrio cholerae infection






## **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ATP6V0D2 (Cat# [RC209007], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ATP6V0D2 antibody (Cat# [TA812166]). Positive lysates [LY403477] (100ug) and [LC403477] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from HEPG2 cell lines by using anti-ATP6V0D2 monoclonal antibody (1:500).





Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-ATP6V0D2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA812166]) (1:2000)



Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-ATP6V0D2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA812166]) (1:2000)