

Product datasheet for **TA812216BM**


MYO18A Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI3F7]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI3F7
Applications:	IHC, WB
Recommended Dilution:	WB 1:2000, IHC 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 600-1000 of human MYO18A (NP_976063) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	myosin XVIII A
Database Link:	NP_976063 Entrez Gene 399687 Human Q92614
Synonyms:	MAJN; MYSPDZ; SP-R210; SPR210



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY MYO18A (Cat# [RC224962], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MYO18A (Cat# [TA812216])(1:2000).



Immunohistochemical staining of paraffin-embedded Human gastric tissue within the normal limits using anti-MYO18A mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA812216]) (1:1000)