

Product datasheet for **TA812221BM**


FUBP3 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI7B6]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI7B6
Applications:	IHC
Recommended Dilution:	IHC 1:500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 2-572 of human FUBP3 (NP_003925) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	61.5 kDa
Gene Name:	far upstream element binding protein 3
Database Link:	NP_003925 Entrez Gene 320267 Mouse Entrez Gene 362106 Rat Entrez Gene 8939 Human Q96I24
Background:	May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression. [UniProtKB/Swiss-Prot Function]
Synonyms:	FBP3



[View online »](#)

Product images:

Immunohistochemical staining of paraffin-embedded Human Gastric Carcinoma using anti-FUBP3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, [TA812221]) (1:500)