

Product datasheet for **TA812224BM**


FUBP3 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI5H6]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI5H6
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 2-572 of human FUBP3 (NP_003925) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	61.5 kDa
Gene Name:	far upstream element binding protein 3
Database Link:	NP_003925 Entrez Gene 320267 Mouse Entrez Gene 362106 Rat Entrez Gene 8939 Human Q96I24
Background:	May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression. [UniProtKB/Swiss-Prot Function]
Synonyms:	FBP3



[View online »](#)

Product images:


HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FUBP3 ([RC215402], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FUBP3. Positive lysates [LY418344] (100ug) and [LC418344] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 3 cell lines by using anti-FUBP3 monoclonal antibody (1:500).