

## Product datasheet for **TA812252BM**


### **GI24 (C10orf54) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI8B8]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI8B8
Applications:	FC
Recommended Dilution:	FLOW 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human C10ORF54 (NP_071436) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	33.7 kDa
Gene Name:	chromosome 10 open reading frame 54
Database Link:	<a href="#">NP_071436</a> <a href="#">Entrez Gene 64115 Human</a> <a href="#">Q9H7M9</a>
Synonyms:	B7-H5; B7H5; GI24; PP2135; SISP1
Protein Families:	Transmembrane



[View online »](#)

**Product images:**

HEK293T cells transfected with either [RC204547] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-C10orf54 antibody ([TA812252]), and then analyzed by flow cytometry (1:100).