

Product datasheet for TA812520

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MRP2 (ABCC2) Mouse Monoclonal Antibody [Clone ID: OTI5C3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5C3
Applications: IHC, WB

Recommended Dilution: WB 1:500, IHC 1:5000

Reactivity: Human, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1256-1545 of human

ABCC2 (NP_000383) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: ATP binding cassette subfamily C member 2

Database Link: NP 000383

Entrez Gene 25303 RatEntrez Gene 1244 Human

Q92887

Background: Mediates hepatobiliary excretion of numerous organic anions. May function as a cellular

cisplatin transporter. [UniProtKB/Swiss-Prot Function]

Synonyms: ABC30; CMOAT; cMRP; DJS; MRP2


Protein Families: Druggable Genome, Transmembrane

Protein Pathways: ABC transporters







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ABCC2 (Cat# [RC218976], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ABCC2 (Cat# TA812520)(1:500).



Western blot analysis of extracts (35ug) from rat liver lysates by using anti-ABCC2 monoclonal antibody (1:500).



Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-ABCC2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA812520) (1:5000)