

## Product datasheet for **TA812559**

### **TIM 3 (HAVCR2) Mouse Monoclonal Antibody [Clone ID: OTI2E5]**

#### **Product data:**

|                         |                                                                                                           |
|-------------------------|-----------------------------------------------------------------------------------------------------------|
| Product Type:           | Primary Antibodies                                                                                        |
| Clone Name:             | OTI2E5                                                                                                    |
| Applications:           | FC, WB                                                                                                    |
| Recommended Dilution:   | WB 1:2000, FLOW 1:100                                                                                     |
| Reactivity:             | Human                                                                                                     |
| Host:                   | Mouse                                                                                                     |
| Isotype:                | IgG2b                                                                                                     |
| Clonality:              | Monoclonal                                                                                                |
| Immunogen:              | Full length human recombinant protein of human HAVCR2 (NP_116171) produced in HEK293T cell.               |
| Formulation:            | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                                      |
| Concentration:          | 1 mg/ml                                                                                                   |
| Purification:           | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Conjugation:            | Unconjugated                                                                                              |
| Storage:                | Store at -20°C as received.                                                                               |
| Stability:              | Stable for 12 months from date of receipt.                                                                |
| Predicted Protein Size: | 33.2 kDa                                                                                                  |
| Gene Name:              | hepatitis A virus cellular receptor 2                                                                     |
| Database Link:          | <a href="#">NP_116171</a><br><a href="#">Entrez Gene 84868 Human</a><br><a href="#">Q8TDQ0</a>            |

[View online »](#)

**Background:**


The protein encoded by this gene belongs to the immunoglobulin superfamily, and TIM family of proteins. CD4-positive T helper lymphocytes can be divided into types 1 (Th1) and 2 (Th2) on the basis of their cytokine secretion patterns. Th1 cells are involved in cell-mediated immunity to intracellular pathogens and delayed-type hypersensitivity reactions, whereas, Th2 cells are involved in the control of extracellular helminthic infections and the promotion of atopic and allergic diseases. This protein is a Th1-specific cell surface protein that regulates macrophage activation, and inhibits Th1-mediated auto- and alloimmune responses, and promotes immunological tolerance. [provided by RefSeq, Sep 2011]

**Synonyms:**


CD366; HAVcr-2; KIM-3; SPTCL; Tim-3; TIM3; TIMD-3; TIMD3

**Protein Families:**



Druggable Genome, Transmembrane

**Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HAVCR2 ([RC209440], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HAVCR. Positive lysates [LY403200] (100ug) and [LC403200] (20ug) can be purchased separately from OriGene.



Flow cytometric analysis of living 293T cells transfected with HAVCR2 overexpression plasmid ([RC209440]), Red/empty vector ([PS100001], Blue) using anti-HAVCR2 antibody (TA812559). Cells incubated with a non-specific antibody (Green) were used as isotype control (1:100).



Flow cytometric analysis of living PBMCs treated with 10ug/ml PHA for 72h (Right)/untreated (Left) using anti-HAVCR2 antibody (TA812559) (1:100).



Flow cytometric analysis of living PBMCs treated with 10ug/ml PHA for 72h (Right) using anti-HAVCR2 antibody (TA812559). Cells incubated with a non-specific antibody (Left) were used as isotype control (1:100).