

## Product datasheet for **TA812736BM**


### ZNF226 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI8G4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI8G4
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ZNF226 (NP_001027545) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	91.7 kDa
Gene Name:	zinc finger protein 226
Database Link:	<a href="#">NP_001027545</a> <a href="#">Entrez Gene 7769 Human</a> <a href="#">Q9NYT6</a>
Background:	May be involved in transcriptional regulation. [UniProtKB/Swiss-Prot Function]
Protein Families:	Transcription Factors



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZNF226 ([RC223775], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF226 (1:500).