

## Product datasheet for **TA812791AM**


### **VGLL3 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI3C7]**

#### **Product data:**


Product Type:	Primary Antibodies
Clone Name:	OTI3C7
Applications:	WB
Recommended Dilution:	WB 1:250~500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human VGLL3 (NP_057290) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	36 kDa
Gene Name:	vestigial like family member 3
Database Link:	<a href="#">NP_057290</a> <a href="#">Entrez Gene 73569 Mouse</a> <a href="#">Entrez Gene 498038 Rat</a> <a href="#">Entrez Gene 389136 Human</a> <a href="#">A8MV65</a>
Synonyms:	VGL-3; VGL3
Protein Families:	Transcription Factors



[View online »](#)

**Product images:**


HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY VGLL3 (Cat# [RC218896], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-VGLL3 (Cat# [TA812791])(1:500).



Western blot analysis of extracts (35ug) from 8 different cell lines by using anti-VGLL3 monoclonal antibody (1:250).