

## Product datasheet for **TA813058BM**

### ITGB1BP3 (NMRK2) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI3H9]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3H9
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ITGB1BP3 (NP_733778) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	nicotinamide riboside kinase 2
Database Link:	<a href="#">NP_733778</a> <a href="#">Entrez Gene 27231 Human</a> <a href="#">Q9NPI5</a>
Synonyms:	ITGB1BP3; MIBP; NRK2



[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ITGB1BP3 (Cat# [RC210779], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ITGB1BP3 (Cat# [TA813058])(1:2000).