

### **Product datasheet for TA813207**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CD70 Mouse Monoclonal Antibody [Clone ID: OTI13D7]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI13D7

Applications: FC

Recommended Dilution: FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human CD70 (NP\_001243) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 20.9 kDa

**Gene Name:** CD70 molecule

Database Link: NP 001243

Entrez Gene 970 Human

P32970

**Background:** Cytokine that binds to CD27. Plays a role in T-cell activation. Induces the proliferation of

costimulated T-cells and enhances the generation of cytolytic T-cells. [UniProtKB/Swiss-Prot

Function]

Synonyms: CD27-L; CD27L; CD27LG; LPFS3; TNFSF7; TNLG8A


**Protein Families:** ES Cell Differentiation/IPS, Transmembrane






**Protein Pathways:** Cytokine-cytokine receptor interaction

# **Product images:**




Flow cytometric analysis of living 293T cells transfected with CD70 overexpression plasmid ([RC200410]), Red)/empty vector ([PS100001], Blue) using anti-CD70 antibody (TA813207). Cells incubated with a non-specific antibody (Green) were used as isotype control (1:100).





Flow cytometric analysis of living Raji cells, using anti-CD70 antibody (TA813207, Red), compared to an isotype control (green), and a PBS control (blue) (1:100).







Flow cytometric analysis of living 786-O cells, using anti-CD70 antibody (TA813207, Red), compared to an isotype control (green), and a PBS control (blue) (1:100).





Flow cytometric analysis of living PBMCs treated with 10ug/ml PHA for 72h (Right)/untreated (Left) using anti-CD70 antibody (TA813207) (1:100).





Flow cytometric analysis of living PBMCs treated with 10ug/ml PHA for 72h using anti-CD70 antibody (TA813207) (Right). Cells incubated with a non-specific antibody (Left) were used as isotype control (1:100).