

Product datasheet for **TA813476S**

CD9 Mouse Monoclonal Antibody [Clone ID: OTI3D5]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3D5
Applications:	FC
Recommended Dilution:	FLOW 1:100
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CD9 (NP_001760) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	25.4 kDa
Gene Name:	CD9 molecule
Database Link:	NP_001760 Entrez Gene 12527 Mouse Entrez Gene 24936 Rat Entrez Gene 928 Human P21926



[View online »](#)

Background:

This gene encodes a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Tetraspanins are cell surface glycoproteins with four transmembrane domains that form multimeric complexes with other cell surface proteins. The encoded protein functions in many cellular processes including differentiation, adhesion, and signal transduction, and expression of this gene plays a critical role in the suppression of cancer cell motility and metastasis. [provided by RefSeq, Jan 2011]

Synonyms:


BTCC-1; DRAP-27; MIC3; MRP-1; TSPAN-29; TSPAN29

Protein Families:


Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Transmembrane

Protein Pathways:


Hematopoietic cell lineage

Product images:


Flow cytometric analysis of living HeLa cells, using anti-CD9 antibody([TA813476], Red), compared to an isotype control (green), and a PBS control (blue).(1:100)




Flow cytometric analysis of living K562 cells, using anti-CD9 antibody([TA813476], Red), compared to an isotype control (green), and a PBS control (blue).(1:100)



Flow cytometric analysis of living RPMI-8226 cells, using anti-CD9 antibody([TA813476], Red), compared to an isotype control (green), and a PBS control (blue).(1:100)



Flow cytometric analysis of living HUVEC cells, using anti-CD9 antibody([TA813476], Red), compared to an isotype control (green), and a PBS control (blue).(1:100)



Flow cytometric analysis of living U87MG cells, using anti-CD9 antibody([TA813476], Red), compared to an isotype control (green), and a PBS control (blue).(1:100)