

## Product datasheet for **TA813739**

### MYBPC3 Mouse Monoclonal Antibody [Clone ID: OTI7H3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI7H3
Applications:	ELISA
Recommended Dilution:	ELISA 1:5000-10000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human MYBPC3 C0-C1 (NP_000247) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Predicted Protein Size:	140.8 kDa
Gene Name:	myosin binding protein C, cardiac
Database Link:	<a href="#">NP_000247</a> <a href="#">Entrez Gene 4607 Human</a> <a href="#">Q14896</a>
Background:	MYBPC3 encodes the cardiac isoform of myosin-binding protein C. Myosin-binding protein C is a myosin-associated protein found in the cross-bridge-bearing zone (C region) of A bands in striated muscle. MYBPC3, the cardiac isoform, is expressed exclusively in heart muscle. Regulatory phosphorylation of the cardiac isoform in vivo by cAMP-dependent protein kinase (PKA) upon adrenergic stimulation may be linked to modulation of cardiac contraction. Mutations in MYBPC3 are one cause of familial hypertrophic cardiomyopathy. [provided by RefSeq, Jul 2008]
Synonyms:	CMD1MM; CMH4; cMyBP-C; FHC; LVNC10; MYBP-C



[View online »](#)

**Protein Families:** Druggable Genome

**Protein Pathways:** Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM)